

Clean Dose Presoak Tablet

Enzymatic formula breaks down protein and starch films

Multi-use and effective

- Presoaking flatware results in consistently clean and spotless results with every wash cycle
- Saves time and labor by helping to reduce need for rewash

User-friendly design

- Safer product handling
- Requires no PPE

Easy on the environment

- Concentrated solid tablets contribute to reduced water, energy, and plastic waste

DIRECTIONS FOR USE

1. Tear open a packet and drop the tablet into a presoak tub. Wash hands after handling product.
2. Fill up the tub with warm water (110° – 120° F) and allow the tablet to dissolve.
3. Place dirty flatware into the solution and allow to soak.



The tablet's solid chemistry is formulated to save on packaging.

99.9%
reduction in plastic waste¹

¹ Reduction in plastic waste is an approximation assuming use of one case of product per month over 12 months, compared with the equivalent use of Ecolab liquid products, which is enabled by the smaller, lighter package size for solid tablets.

CLEAN DOSE PRESOAK TABLET

PACK SIZE	48 ct
DILUTION RATIO	1 tablet/tub
MFR #	6102788
SHAMROCK #	

We thought small so you can think big.

Ready to learn more about our full tablet cleaning program?
Visit procleansouthwest.com/cleandose or scan the QR code.

The experienced and dedicated ProClean sales team provides you personalized support to help you run a clean and safe operation.

Call **1 800 35 CLEAN** or visit procleansouthwest.com for more information.

