PROCLEAN

TB Disinfectant Cleaner Ready-To-Use

EPA Reg. No. 1839-83-1677

EPA-registered cleaner disinfects and deodorizes hard, non-porous surfaces with broad-spectrum efficacy.

- Kills SARS-CoV-2 in 60 seconds¹
- Kills Norovirus in 30 second¹
- Convenient, ready-to-use cleaner, disinfectant and deodorizer
- Use on hard, non-porous surfaces such as porcelain ceramic tile, metal, stainless steel and plastic.
- Citrus scent leaves restrooms and guest areas smelling clean and fresh

Directions: (see label for complete directions)

- 1. Remove visible soil prior to application of the product.
- 2. Hold container 6 to 8 inches from surface.
- 3. Spray area until it is covered with the solution.
- 4. Allow product to penetrate and remain wet for 3 minutes.²
- 5. Wipe off with a clean cloth, mop or sponge.

1 When used as directed.

 $2\ {\rm Allow}\ {\rm surface}\ {\rm to}\ {\rm remain}\ {\rm wet}\ {\rm for}\ 5\ {\rm minutes}\ {\rm for}\ {\rm HBV}\ {\rm and}\ {\rm HCV}.$

Housekeeping Solutions

PRODUCT	PACK SIZE	MFR #	SHAMROCK #
TB Disinfectant Cleaner Ready-To-Use	6 - 32 oz	6101593	4136931





The ProClean experienced and dedicated sales team provide you personalized support to help you run a clean and safe operation.

Call 1 800 35 CLEAN or visit procleansouthwest.com for more information.



